# **MPC-5V51**5**D**

### 2~8 ℃ Pharmacy Refrigerator





## **FEATURERS**



### Max. safety for storage

Dual cooling system, dual assurance for



### **Footprint unlock**

More convenient and safer



### **Universal castor wheels**

Proven heavy-duty lockable castors



### Suction circulation system

Less temp. flutuation & better uniformity




### Fast pull down

Real forced air cooling system



### 2 temp. simulation boxes

Accurately detect the temperature of



Mether Dural CoolingSystem



### 3 levels of authority

Multiple application scenarios



#### **NFC** unlock

Support mobile phone NFC and IC card



### **Full information display**

Alarms, curve, temperature and permission



### 8 types of alarm

Audible and visual alarms



### Precise temp. control

Temp. control accurate to 0.1°C



### **Eco-friendly**

CFC-free refrigerant and energy-saving



# **STRUCTURE**





High/low temp. alarm

Condenser fault alarm

System abnormal alarm

Thermostat failure alarm

Door ajar alarm

Power failure alarm

**Ambient Temperature** 





Real forced air circulation and excellent refrigeration system pull down the chamber temperature rapidly. And keep the chamber with good temperature uniformity and fluctuation.

# **PARAMETER**

Performance		Refrigeration System		Alarms	
Cooling method	Forced Air Cooling	Refrigerant	R600a	High/low temperature	Υ
Defrost	Auto	Temp. controller	Microprocessor	Power failure	Υ
Shelf/drawer	2/4	Sensor	NTC	Thermostat failure	Υ
Temperature range	2°C~8°C	Display	Digital display	System abnormal alarm	Υ
Ambient temperature	10~32°C	External material	SPCC sprayed steel	Door ajar	Υ
NT./GT.	140/157	Capacity (L/cu.ft)	515	Ambient Temperature	Υ
Noise (db)	50	Compressor QTY	EMBRACO/2	Condenser fault	Υ
Electrical Data		Dimensions		Other	
Voltage frequency( V/Hz)	220V 50Hz/60Hz	Package size(W*D*H)(mm)	755*855*2115	Inside door	2 shelves + 4 drawers
Power (W)	294W	Exterior size(W*D*H)(mm)	695*794*1993	Remote alarm port	Optional
Current (A)	2.1A	Interior size(W*D*H)(mm)	595*641*1475	Access port QTY/Dia.	1/25mm
Power consumption	3.9kWh/24h	СВМ	1.37	Internal material	SPCC sprayed steel